

Product Datasheet

C1QTNF6 antibody (FITC) (orb37429)

Description	C1QTNF6 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Complement C1q tumor necrosis factor-related protein 6 protein (47-278AA)
Target	C1QTNF6
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9BXI9
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.