


Product Datasheet

GABRA5 antibody (orb374299)

Description	Rabbit polyclonal antibody to GABRA5
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human GABRA5(NP_000801.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	52kDa
Uniprot ID	P31644
Dilution Range	WB: 1:1000-4000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...