


Product Datasheet

GLDC antibody (orb374300)

Description	Rabbit polyclonal antibody to GLDC
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, IP, WB
Immunogen	Recombinant fusion protein of human GLDC(NP_000161.2).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	113kDa
Uniprot ID	P23378
Dilution Range	WB: 1:200-500
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...