


Product Datasheet

SAE1 antibody (orb374315)

Description	Rabbit polyclonal antibody to SAE1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human SAE1 (NP_005491.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	38kDa
Uniprot ID	Q9UBE0
Entrez	10055
Dilution Range	WB: 1:200-500
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...