

Product Datasheet

KRTAP3-3 antibody (orb37432)

Description	KRTAP3-3 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human Keratin-associated protein 3-3 protein (1-98AA)
Target	KRTAP3-3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9BYR6
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.