


Product Datasheet

PABPC5 antibody (orb374336)

Description	Rabbit polyclonal antibody to PABPC5
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human PABPC5(NP_543022.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	43kDa
Uniprot ID	Q96DU9
Dilution Range	WB: 1:1000-4000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...