


Product Datasheet

FYN (Phospho-S21) antibody (orb374340)

Description	Rabbit polyclonal antibody to FYN (Phospho-S21)
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around S21 of human Phospho-Fyn-S21(NP_002028.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	60kDa
Uniprot ID	P06241
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of A47...