


Product Datasheet

GSK3A (Phospho-S21) antibody (orb374350)

Description	Rabbit polyclonal antibody to GSK3A (Phospho-S21)
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	A synthetic phosphorylated peptide around S21 of human GSK3 α (NP_063937.2).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	51kDa
Uniprot ID	P49840
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...