


Product Datasheet

CHEK2 (Phospho-T68) antibody (orb374351)

Description	Rabbit polyclonal antibody to CHEK2 (Phospho-T68)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around T68 of human Chk2 (NP_009125.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	62kDa
Uniprot ID	O96017
Entrez	11200
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of 293...