


Product Datasheet

BCR (Phospho-Y177) antibody (orb374353)

Description	Rabbit polyclonal antibody to BCR (Phospho-Y177)
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around Y177 of human BCR(NP_004318.3).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	160kDa
Uniprot ID	P11274
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of K562...