


Product Datasheet

CYLD antibody (orb374370)

Description	Rabbit polyclonal antibody to CYLD
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human CYLD(NP_001035877.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	107KDa
Uniprot ID	Q9NQC7
Dilution Range	WB: 1:200-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...