


Product Datasheet


UBL3 antibody (orb374374)

Description	Rabbit polyclonal antibody to UBL3
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human UBL3(NP_009037.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	13kDa
Uniprot ID	O95164
Dilution Range	WB: 1:1000-2000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of C6 cells

...



Western blot analysis of extracts of var...