



Product Datasheet

DHX38 antibody (orb374382)

Description	Rabbit polyclonal antibody to DHX38
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human DHX38(NP_054722.2).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	140kDa
Uniprot ID	Q92620
Dilution Range	WB: 1:1000-2000, IP: 1:20-50
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human gastric tissue.



Western blot analysis of extracts of various cell lines.