


Product Datasheet

OSBPL3 antibody (orb374390)

Description	Rabbit polyclonal antibody to OSBPL3
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human OSBPL3(NP_056365.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	101kDa
Uniprot ID	Q9H4L5
Dilution Range	WB: 1:200-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...