


Product Datasheet

BLOC1S6 antibody (orb374391)

Description	Rabbit polyclonal antibody to BLOC1S6
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human BLOC1S6(NP_036520.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	20kDa
Uniprot ID	Q9UL45
Dilution Range	WB: 1:200-2000, IHC-P: 1:20-200, IP: 1:20-50
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...