



Product Datasheet

SUZ12 antibody (orb374425)


Description	Rabbit polyclonal antibody to SUZ12
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IP, WB
Immunogen	Recombinant fusion protein of human SUZ12(NP_056170.2).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	83kDa
Uniprot ID	Q15022
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rat hear...



Immunohistochemical staining of rat brai...



Immunohistochemical staining of rat lung...