



Product Datasheet

GDI2 antibody (orb374455)


Description	Rabbit polyclonal antibody to GDI2
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human GDI2(NP_001485.2).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	51kDa
Uniprot ID	P50395
Dilution Range	WB: 1:200-3000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of NIH/3T3 C...



Western blot analysis of extracts from n...



Western blot analysis of extracts of var...