


Product Datasheet

VPS33B antibody (orb374470)

Description	Rabbit polyclonal antibody to VPS33B
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human VPS33B(NP_061138.3).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	62kDa
Uniprot ID	Q9H267
Dilution Range	WB: 1:1000-2000, IHC-P: 1:20-200, IF/ICC: 1:10-100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...