


Product Datasheet

SIRPA antibody (orb374479)

Description	Rabbit polyclonal antibody to SIRPA
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human CD172a/SIRPα(NP_542970.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	65kDa
Uniprot ID	P78324
Dilution Range	WB: 1:200-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...