



## Product Datasheet

### ZFYVE9 antibody (orb374481)

|                            |                                                              |
|----------------------------|--------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to ZFYVE9                         |
| <b>Species/Host</b>        | Rabbit                                                       |
| <b>Reactivity</b>          | Human, Mouse, Rat                                            |
| <b>Conjugation</b>         | Unconjugated                                                 |
| <b>Tested Applications</b> | WB                                                           |
| <b>Immunogen</b>           | Recombinant fusion protein of human ZFYVE9 (NP_004790.2).    |
| <b>Preservatives</b>       | Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.    |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.                  |
| <b>Note</b>                | For research use only                                        |
| <b>Isotype</b>             | IgG                                                          |
| <b>Clonality</b>           | Polyclonal                                                   |
| <b>Purity</b>              | Affinity purification                                        |
| <b>MW</b>                  | 156kDa                                                       |
| <b>Uniprot ID</b>          | <b>O95405</b>                                                |
| <b>Entrez</b>              | <b>9372</b>                                                  |
| <b>Dilution Range</b>      | WB: 1:200-2000, IF/ICC: 1:10-100, IP: 1:20-50, FACS: 1:20-50 |
| <b>Expiration Date</b>     | 12 months from date of receipt.                              |



Immunofluorescence analysis of HeLa cell...



Western blot analysis of extracts of var...