

Product Datasheet

IMPDH1 antibody (orb374496)

Description	Rabbit polyclonal antibody to IMPDH1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human IMPDH1(NP_899066.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	55KDa
Uniprot ID	P20839
Dilution Range	WB: 1:1000-2000
Expiration Date	12 months from date of receipt.