



Product Datasheet

UBQLN2 antibody (orb374500)


Description	Rabbit polyclonal antibody to UBQLN2
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human UBQLN2(NP_038472.2).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	70kDa
Uniprot ID	Q9UHD9
Dilution Range	WB: 1:200-2000, IHC-P: 1:20-200, IF/ICC: 1:20-100
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of NIH-3T3 C...



Immunofluorescence analysis of HeLa cell...



Western blot analysis of extracts of var...