


Product Datasheet

UACA antibody (orb374503)

Description	Rabbit polyclonal antibody to UACA
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human UACA (NP_060473.2).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	162kDa
Uniprot ID	Q9BZF9
Entrez	55075
Dilution Range	WB: 1:500-2000, IF/ICC: 1:20-100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...