



Product Datasheet

PMAIP1 antibody (orb374505)


Description	Rabbit polyclonal antibody to PMAIP1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	A synthetic peptide of human PMAIP1(NP_066950.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	14kDa
Uniprot ID	Q13794
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of mouse bre...



Immunofluorescence analysis of mouse bre...



Immunofluorescence analysis of rat breas...