


Product Datasheet

Syk (Phospho-Y525/526) antibody (orb374525)

Description	Rabbit polyclonal antibody to Syk (Phospho-Y525/526)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around Y525 & Y526 of human Syk (NP_003168.2).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	72kDa
Uniprot ID	P43405
Entrez	6850
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of Jur...