

Product Datasheet

ACK1 (Phospho-Y284) antibody (orb374529)

Description	Rabbit polyclonal antibody to ACK1 (Phospho-Y284)
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around Y284 of human TNK2(NP_005772.3).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	180kDa
Uniprot ID	Q07912
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.