


Product Datasheet

Golgin subfamily A member 4 antibody (orb374539)

Description	Rabbit polyclonal antibody to Golgin subfamily A member 4.
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human GOLGA4(NP_002069.2).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	280kDa
Uniprot ID	Q13439
Dilution Range	WB: 1:1000-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...