



Product Datasheet

HINT1 antibody (orb374542)

Description	Rabbit polyclonal antibody to HINT1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human HINT1(NP_005331.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	14kDa
Uniprot ID	P49773
Dilution Range	WB: 1:1000-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from n...



Western blot analysis of extracts of var...