


Product Datasheet

SLC11A2 antibody (orb374547)

Description	Rabbit polyclonal antibody to SLC11A2
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human DMT1/SLC11A2(NP_001167597.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	72KDa
Uniprot ID	P49281
Dilution Range	WB: 1:1000-2000
Expiration Date	12 months from date of receipt.



Western blot
analysis of
extracts of
var...