


Product Datasheet

XPNPEP2 antibody (orb374556)

Description	Rabbit polyclonal antibody to XPNPEP2
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human XPNPEP2 (NP_003390.4).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	76kDa
Uniprot ID	043895
Entrez	7512
Dilution Range	WB: 1:1000-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...