



Product Datasheet

KCNMB2 antibody (orb374569)

Description	Rabbit polyclonal antibody to KCNMB2
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human KCNMB2(NP_005823.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	30kDa
Uniprot ID	Q9Y691
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human apical surface of epithelial cells.



Western blot analysis of extracts of various tissues.