

Product Datasheet

AMHR2 antibody (orb37457)

Description	Rabbit polyclonal antibody to AMHR2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	AMHR2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	62750
Uniprot ID	Q16671
NCBI	NP_001158162.1 , NP_065434.1 , NP_001158163.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

