




Product Datasheet

AMHR2 antibody (orb37458)

Description	Rabbit polyclonal antibody to AMHR2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Synthetic Peptide
Target	AMHR2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	62750
Uniprot ID	Q16671
NCBI	NP_001158162.1 , NP_065434.1 , NP_001158163.1
Dilution Range	WB: 1:1000, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Western blot analysis of 293 cell lysate...



Flow cytometric analysis of CEM cells us...