


Product Datasheet

CABLES1 antibody (orb374594)

Description	Rabbit polyclonal antibody to CABLES1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human CABLES1(NP_612384.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	67kDa
Uniprot ID	Q8TDN4
Dilution Range	WB: 1:1000-2000
Expiration Date	12 months from date of receipt.



Western blot
analysis of
extracts of
var...