


Product Datasheet

GORASP2 antibody (orb374612)

Description	Rabbit polyclonal antibody to GORASP2
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant fusion protein of human GORASP2 (NP_056345.3).
Preservatives	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	55kDa
Uniprot ID	Q9H8Y8
Entrez	26003
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.



Western blot
analysis of
extracts of
var...