



Product Datasheet

UBE2G2 antibody (orb374622)

Description	Rabbit polyclonal antibody to UBE2G2
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human UBE2G2(NP_003334.2).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	19kDa
Uniprot ID	P60604
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of NIH-3T3 C...



Western blot analysis of extracts of var...