



Product Datasheet

NFKBIZ antibody (orb374656)

Description	Rabbit polyclonal antibody to NFKBIZ
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human NFKBIZ(NP_113607.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	64KDa
Uniprot ID	Q9BYH8
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of Raw264.7 ...



Western blot analysis of extracts of var...