


Product Datasheet

SIGLEC8 antibody (orb374665)

Description	Rabbit polyclonal antibody to SIGLEC8
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human Siglec-8(NP_055257.2).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	60kDa
Uniprot ID	Q9NYZ4
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...