



Product Datasheet

NOTCH2 antibody (orb374673)


Description	Rabbit polyclonal antibody to NOTCH2
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human NOTCH2(NP_077719.2).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	120kDa
Uniprot ID	Q04721
Dilution Range	WB: 1:500-2000, IHC-P: 1:100-200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...



Immunohistochemical staining of human st...



Immunohistochemical staining of mouse br...