


Product Datasheet

Collagen III antibody (orb374676)

Description	Rabbit polyclonal antibody to Collagen III
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human III alpha 1/COL3A1(NP_000081.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	200KDa
Uniprot ID	P02461
Dilution Range	WB: 1:500-2000, IHC-P: 1:20-200, IF/ICC: 1:20-100, FACS: 1:20-50
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...