


## Product Datasheet

### MRP 2 antibody (orb374693)

|                            |                                                              |
|----------------------------|--------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to MRP 2                          |
| <b>Reactivity</b>          | Human, Mouse, Rat                                            |
| <b>Conjugation</b>         | Unconjugated                                                 |
| <b>Tested Applications</b> | ICC, IF, WB                                                  |
| <b>Immunogen</b>           | Recombinant fusion protein of human MRP2/ABCC2(NP_000383.2). |
| <b>Preservatives</b>       | PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.             |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.                  |
| <b>Note</b>                | For research use only                                        |
| <b>Isotype</b>             | IgG                                                          |
| <b>Clonality</b>           | Polyclonal                                                   |
| <b>MW</b>                  | 200KDa-270KDa                                                |
| <b>Uniprot ID</b>          | <b>Q92887</b>                                                |
| <b>Dilution Range</b>      | WB: 1:500-2000                                               |
| <b>Expiration Date</b>     | 12 months from date of receipt.                              |



Western blot analysis of extracts of var...