

Product Datasheet

CD105 peptide (orb374706)

Description	CD105 peptide
Conjugation	Unconjugated
Tested Applications	HPLC
Target	CD105
Form/Appearance	Lyophilized
Storage	Upon reconstitution, aliquot and store at -20°C or -80°C. for long term storage, Avoid repeated freeze/thaw cycles.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb10285 (CD105 antibody)
Source	Synthesis
Expiration Date	6 months from date of receipt.