

Product Datasheet

Melatonin Receptor 1A Peptide (orb374736)

Description	Melatonin Receptor 1A Peptide
Conjugation	Unconjugated
Tested Applications	HPLC
Target	Melatonin Receptor 1A
Form/Appearance	Lyophilized
Storage	Upon reconstitution, aliquot and store at -20°C or -80°C. for long term storage, Avoid repeated freeze/thaw cycles.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb11085 (Melatonin Receptor 1A antibody)
Source	Synthesis
Expiration Date	6 months from date of receipt.