

Product Datasheet

CD163 peptide (orb374759)

Description	CD163 peptide
Conjugation	Unconjugated
Tested Applications	HPLC
Target	CD163
Form/Appearance	Lyophilized
Storage	Upon reconstitution, aliquot and store at -20°C or -80°C. for long term storage, Avoid repeated freeze/thaw cycles.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb13303 (CD163 antibody)
Source	Synthesis
Expiration Date	6 months from date of receipt.