

Product Datasheet

Ly-6G peptide (orb374770)

Description	Ly-6G peptide
Conjugation	Unconjugated
Tested Applications	HPLC
Target	Ly-6G
Form/Appearance	Lyophilized
Storage	Upon reconstitution, aliquot and store at -20°C or -80°C. for long term storage, Avoid repeated freeze/thaw cycles.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb13552 (Ly-6G antibody)
Source	Synthesis
Expiration Date	6 months from date of receipt.