

Product Datasheet

PGC1 alpha peptide (orb374773)

Description	PGC1 alpha peptide
Conjugation	Unconjugated
Tested Applications	HPLC
Target	PGC1 alpha
Form/Appearance	Lyophilized
Storage	Upon reconstitution, aliquot and store at -20°C or -80°C. for long term storage, Avoid repeated freeze/thaw cycles.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb13647 (PGC1 alpha antibody)
Source	Synthesis
Expiration Date	6 months from date of receipt.