

Product Datasheet

Prostaglandin E synthase 2 peptide (orb374831)

Description	Prostaglandin E synthase 2 peptide
Conjugation	Unconjugated
Tested Applications	HPLC
Target	Prostaglandin E synthase 2
Form/Appearance	Lyophilized
Storage	Upon reconstitution, aliquot and store at -20°C or -80°C. for long term storage, Avoid repeated freeze/thaw cycles.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb191502 (Prostaglandin E synthase 2 antibody)
Source	Synthesis
Expiration Date	6 months from date of receipt.