

Product Datasheet

Integrin beta 1 peptide (orb374836)

Description	Integrin beta 1 peptide
Conjugation	Unconjugated
Tested Applications	HPLC
Target	Integrin beta 1
Form/Appearance	Lyophilized
Storage	Upon reconstitution, aliquot and store at -20°C or -80°C. for long term storage, Avoid repeated freeze/thaw cycles.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb197998 (Integrin beta 1 antibody)
Source	Synthesis
Expiration Date	6 months from date of receipt.