

## Product Datasheet

### beta 1 Adrenergic Receptor peptide (orb374842)

|                            |                                                                                                                     |
|----------------------------|---------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | beta 1 Adrenergic Receptor peptide                                                                                  |
| <b>Conjugation</b>         | Unconjugated                                                                                                        |
| <b>Tested Applications</b> | HPLC                                                                                                                |
| <b>Target</b>              | beta 1 Adrenergic Receptor                                                                                          |
| <b>Form/Appearance</b>     | Lyophilized                                                                                                         |
| <b>Storage</b>             | Upon reconstitution, aliquot and store at -20°C or -80°C. for long term storage, Avoid repeated freeze/thaw cycles. |
| <b>Note</b>                | For research use only                                                                                               |
| <b>Application notes</b>   | This peptide was used as immunogen to produce orb229828 (beta 1 Adrenergic Receptor antibody)                       |
| <b>Source</b>              | Synthesis                                                                                                           |
| <b>Expiration Date</b>     | 6 months from date of receipt.                                                                                      |