

Product Datasheet

Integrin alpha 4 peptide (orb374845)

Description	Integrin alpha 4 peptide
Conjugation	Unconjugated
Tested Applications	HPLC
Target	Integrin alpha 4
Form/Appearance	Lyophilized
Storage	Upon reconstitution, aliquot and store at -20°C or -80°C. for long term storage, Avoid repeated freeze/thaw cycles.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb229830 (Integrin alpha 4 antibody)
Source	Synthesis
Expiration Date	6 months from date of receipt.